



## Process Mass Spectrometer, Model AIRSENSE

### On-Line Process Gas Analysis

#### Fast Response and Selectivity

Stricter and stricter environmental laws and regulations has resulted in new and elaborate technologies.

Knowing the customers high profile expectations combined with more than 15 years of experience, the model **AIRSENSE** reflects the latest developmental step within the series of our process mass spectrometer.

Focusing applications for non limited gas components, on line gas analysis should offer reliability, highest possible flexibility, and robustness all in once.

The **AIRSENSE** is based on the Ion Molecule Reaction (IMR) mass spectrometer principle. Unlike others, the IMR mass spectroscopy offers an unique fast response time, a wide dynamic range, selective measurements and lowest detectable limits.


#### Applications Range widely

Based on the principle of measurement and the unique application capability ranges from automotive and catalytic developments, fuel cell R&D, waste burners and power plants, the food industry, pharmaceutical R&D and medical diagnostics.

#### Proven Technology

Minimized service and operational costs, flexibility and user friendly interface are achieved by the combining well proven technology and state of the art manufacturing capabilities.

A user friendly software package contains system configurations and data reporting issues being driven by either by an internal PC-Controller or an external AK- or PLC master.



- |                               |                             |                     |
|-------------------------------|-----------------------------|---------------------|
| 1. Primary Ion Source         | 3. Charge Exchange Cell     | 6. Vacuum System    |
| 2. Octopole Separation Device | 4. Quadrupole - Mass Filter | 7. Gas Inlet System |
| 5. Particle Detector          |                             |                     |

## Principle of Operation

Targeting the specific needs within the wide range of applications, the series AIRSENSE is based on the measurement principle called Ion Molecule Reaction (IMR).

By using the IMR technique, the measurement becomes fast, selective and free of any interference. Unlike others, no fragmentation or overlapping spectra can damage the interpretation of the detected results.

IMR means using primary ions with lower energy level (10 eV .. 14 eV) to completely ionize the probe gas molecules. The signal/noise ratio will be optimized by the integrated octopole separator, focusing the primary ions and filtering out any interferences. The quadrupole mass filter (0-500 amu) then separates the molecules for further detection at the fast pulse counter.

The temperature and pressure compensated sample gas inlet guarantees the correct measurement and avoids any discrimination of the gas.

Any contamination due to condensation or particulate matters is minimized.

## Features Functions Benefits

- Multi component analyzer, high selectivity using IMR technology
  - Wide dynamic range, lowest LDL, linear across the whole range setting
  - Fast response time, optimized operational cost

## Technical details

Mass range	0 – 500 amu	Ambient temperature	20°C - 40°C,
Resolution	< 1 amu	Temperature change	max.1°C/h
Analysis time	10 – 6500 msec/amu	Humidity	max 80% (none condensing)
Measuring range	10 <sup>4</sup>	Inlet Temperature	80°C – 190°C adjustable
Response time	T90 < 50 msec	Gas consumption	30 – 250 ml/min
Lower detection limit	< 1 ppb (Benzene in air) < 10 ppb (Benzene in exhaust gas)	Power supply	220V/50Hz (115V/60Hz), 800 W
Drift concentration	< ± 5% over 12 h (1ppm Benzene)	Dimensions (WxHxD)	590 x 650 x 730 mm
Reproducibility	< ± 3% (1ppm Benzene)	Weight	87 kg
Accuracy	< ± 2% (1ppm Benzene)		

